


Performance of TEM-PCR versus Culture for Bacterial Identification in Pediatric Musculoskeletal Infections

James Wood, MD, MSCI¹, Cheryl Sesler, MS², Donald Stalons, PhD, D(ABMM), MPH², Elena Grigorenko, PhD², C. Buddy Creech, MD, MPH³, Isaac Thomsen, MD, MSCI³

¹Pediatric Infectious Diseases, Indiana University School of Medicine, Indianapolis, IN; ²Diatherix-Eurofins, Huntsville, AL; ³Pediatric Infectious Diseases, Vanderbilt University School of Medicine, Nashville, TN



INTRODUCTION

- Musculoskeletal infections (MSKI) in children require prompt diagnosis and treatment due to risk of local tissue damage and metastatic bacterial spread.
- Despite blood and tissue sampling, cultures remain negative in approximately half of all cases.
- There is a need for improved rapid diagnostics in children with MSKI.
- We compared the detection of pathogens by culture vs. target-enriched multiplex PCR (TEM-PCR™) in children with acute MSKI.

METHODS

- Children 6 months to 18 years of age admitted to Vanderbilt Children's Hospital with MSKI were prospectively enrolled.
- Synovial fluid and bone samples were collected.
- Bacterial cultures and antibiotic susceptibility testing were performed by the Vanderbilt University Medical Center Clinical Laboratory.
- Samples were evaluated by TEM-PCR (Fig. 1) for detection of common pathogens:
 - Staphylococcus aureus* [including methicillin (*mecA*) and clindamycin/erythromycin (*ermA* or *ermC*) resistance genes and the Pantone-Valentine leukocidin (PVL) locus]
 - Kingella kingae*
 - Haemophilus influenzae*
 - Streptococcus pyogenes*
 - Streptococcus pneumoniae*




Figure 1. TEM-PCR™ Scheme
Low concentration nested gene-specific primers (Fo-forward out; Fi-forward in; Ri-reverse in; and Ro-reverse out) are designed to enrich the targets during the initial PCR cycles. A pair of universal SuperPrimers (Fs and Rs) is subsequently used to amplify all targets. The Rs primer is labeled with biotin for amplicon detection.

RESULTS

Table 1. Baseline Characteristics of Patient Population

Characteristic	N=25
Sex	
Male	14 (56)
Female	11 (44)
Race	
Non-Hispanic White	17 (68)
Non-Hispanic Black	6 (24)
Other	2(8)
Age, yr (median [IQR])	8 (3-12)

Data are presented as number (percent) unless otherwise indicated, IQR, interquartile range

RESULTS

Table 2. Comparison of Pathogen Recovery and Antibiotic Resistance Testing For Each Sample

Specimen	Culture Result	Antibiotic Resistance (Culture-based method)			TEM-PCR Result	Antibiotic Resistance Gene Detection (PCR method)			PVL
		Methicillin	Clindamycin	Erythromycin		Methicillin	Clindamycin	Erythromycin	
1	<i>S. aureus</i>	-	+	+	<i>S. aureus</i>	-	+	+	-
2	<i>S. aureus</i>	-	-	-	<i>S. aureus</i>	-	-	-	-
3	No growth				Negative				
4	<i>S. aureus</i>	+	-	-	<i>S. aureus</i>	+	-	-	+
5	<i>S. aureus</i>	+	-	+	<i>S. aureus</i>	+	-	-	+
6	No growth				Negative				
7	<i>S. aureus</i>	-	-	-	<i>S. aureus</i>	-	-	-	+
8	<i>S. aureus</i>	-	-	-	<i>S. aureus</i>	-	-	-	-
9	<i>S. aureus</i>	-	-	-	<i>S. aureus</i>	-	-	-	-
10	<i>S. aureus</i>	-	-	+	<i>S. aureus</i>	-	-	-	-
11	<i>S. aureus</i>	+	-	+	<i>S. aureus</i>	+	-	-	+
12	<i>S. aureus</i>	-	-	-	<i>S. aureus</i>	-	-	-	+
13	<i>S. aureus</i>	-	+	+	<i>S. aureus</i>	-	+	+	-
14	<i>S. aureus</i>	-	-	+	<i>S. aureus</i>	-	-	-	+
15	No growth				Negative				+
16	<i>S. aureus</i>	-	-	-	<i>S. aureus</i>	-	-	-	-
17	<i>S. aureus</i>	-	-	-	<i>S. aureus</i>	-	-	-	-
18	No growth				<i>K. kingae</i>				
19	No growth ^a				<i>K. kingae</i>				
20	No growth ^b				<i>S. aureus</i>	-	-	-	+
21	<i>S. aureus</i>	+	-	+	<i>S. aureus</i>	+	-	-	-
22	<i>S. aureus</i>	+	-	+	<i>S. aureus</i>	+	-	-	+
23	No growth ^c				Negative				
24	No growth ^c				Negative				
25	<i>S. aureus</i>	-	-	-	<i>S. aureus</i>	-	-	-	-

^a Positive for *K. kingae* by PCR, ^b Subject 14 and 20 are the same, samples from different dates, ^c Diagnosed with non-infectious arthritis, PVL- Pantone Valentine Leukocidin




Figure 2. Comparison of Culture vs. TEM-PCR by Specific Pathogen (A) and Antibiotic Resistance (B)




Figure 3. Site of Infection In Children with MSKI

CONCLUSIONS

- TEM-PCR showed 100% concordance with conventional, culture-based testing (gold standard) in pathogen identification.
- TEM-PCR increased pathogen identification to 80%, compared to 68% by culture alone. This effect was seen for both *K. kingae* and *S. aureus*.
- TEM-PCR detection for erythromycin resistance is not inclusive of all known resistance genes, resulting in 11% detection compared to 47% by culture for *S. aureus* specimens.
- TEM-PCR provided rapid antibiotic resistance results that were concordant with culture-based results for methicillin and clindamycin resistance.

IMPLICATIONS

- TEM-PCR has the potential to inform antibiotic selection early in the disease course, decreasing the use of broad-spectrum, empiric antibiotic regimens, and promoting antimicrobial stewardship.
- TEM-PCR, may be a useful adjunct to conventional, culture-based testing for children with acute MSKI.

ACKNOWLEDGEMENTS

We are thankful for the funding support provided by Diatherix-Eurofins Clinical Diagnostics Division.